

YMC30



- C30 chains
- very lipophilic
- exceptional selectivity pattern
- isomer recognition
- polar carotenes
- polar and nonpolar Xanthophylls
- steroids
- retinols
- fat-soluble vitamins
- LC-MS applications



YMC30	Specification
Particle size / μm	3; 5*
Pore size / nm	proprietary
Surface area / m^2g^{-1}	proprietary
Carbon content / %	proprietary
Recommended pH range	2 - 6

* Preparative grade material with 15 μm particle size is available on request.

General

The separation of geometric and positional isomers is a challenging task in reversed phase chromatography. Subtle molecular differences have to be recognized and resolved by the particular stationary phase. Sander et al. have conclusively shown that polymeric C30 HPLC stationary phases are able to discriminate isomeric structures of long chain molecules [1,2].

Properties

Compared to classical C18 stationary phases, YMC30 is much more hydrophobic. Even when pure organic eluents are applied, many sample solutes are retained. The use of non-aqueous reversed phase mobile phases facilitates 100 % solvent recycling and LC-MS applications.

The YMC30 stationary phase provides sufficient phase thickness to enhance interaction with long chained molecules (fig. 1). Therefore, geometric and positional isomers of conjugated double bonding systems are recognized and resolved by the YMC30 phase.

The resolving power of YMC30 for isomers can be verified by the separation of carotenoids, which has been subject of considerable research efforts in the past. Carotenoids are found in a variety of natural sources including fruits and vegetables. In addition, carotenoids are considered as potential drugs for cancer intervention or prevention. Despite the complexity of carotenoid extracts and the minor shape differences between carotenoid isomers, the separation, identification and quantification of these compounds can be achieved by using YMC30 columns.









Fig. 1: Comparison of the film thickness of C18 and C30 stationary phases with the molecular length of β -carotene (determined with Small Angle Neutron Scattering (SANS)).



Applications

YMC30 columns are successfully used in the food industry, for the analysis of vitamin formulations, in environmental analysis, and for the control of the growth of algae. Other potential applications include the separation of prostaglandins or leucotrienes.

Separation of natural compounds



References

- [1] Sander, L.C. and S.A. Wise; *J. Chromatogr.* 1993, 656, 335-351
[2] Sander, L.C. et al.; *Anal. Chem.* 1994, 66, 1667-1674
[3] Block, G. and L. Langseth, "Antioxidant Vitamins and Disease Prevention", *Food Technology* July 1994

^a Courtesy of L.C. Sander, NIST, Gaithersburg, NC, USA

^b Courtesy of J. Schmid, Institut für Seenforschung, Langenargen, Germany



北京慧德易科技有限责任公司

<http://www.herbs-extract.com>

北京慧德易科技有限责任公司致力于为广大提取物厂家、制药企业、科研机构、大专院校以及化妆品、食品、化工企业提供各种天然提取物和一流的实验室仪器、设备以及试剂、耗材等产品。重点服务于广大的天然产物提取物、药物、生物制品的生产厂商，为客户提供从实验室到工业化生产全过程所需的分离设备和耗材，解决您在分离纯化的生产放大的过程中遇到的问题。

我公司代理的产品有罗门哈斯 XAD 系列大孔吸附树脂、三菱化学 MCI Gel, HP、SP 系列树脂产品、Amersham 公司 Sephadex LH-20、MERCK 分析试剂和色谱产品、TOSOH TSK-Gel 和 Toyopearl 凝胶产品、日本 Develosil 和 YMC 高压液相色谱柱和填料，Sigma-Aldrich 标准试剂，以及 Varian、Novasep, Applexsion 等世界著名公司的制备色谱产品，我们除了负责以上各公司产品在中国的销售和市场开发以外，还结合本公司的技术专长，人才优势，积极开展分离技术、色谱技术、样品前处理技术的服务工作。

地址：北京回龙观龙冠大厦 719 室
电话：010-51528296/97/98
Email: sales@herbs-extract.com

邮编：102208
传真：010-51528299
<http://www.herbs-extract.com>